

CLINIPET⁺ FIXED VOLUME PIPETTE SERIES

MODEL	CAT. NO.	COLOUR CODE	VOLUME (μL)	A (%)	P (%)	FIT TO STANDARD TIPS
CP2*	2202-CP	YELLOW	2	± 3.0	± 1.5	200 μL
CP5	2204-CP		5	± 2.0	± 1.5	
CP7*	2206-CP		7	± 1.5	± 0.7	
CP10	2208-CP		10	± 1.2	± 0.5	
CP15*	2210-CP		15	± 1.0	± 0.5	
CP20	2212-CP		20	± 0.9	± 0.4	
CP25	2214-CP		25	± 0.9	± 0.4	
CP30*	2216-CP		30	± 0.9	± 0.3	
CP40*	2218-CP		40	± 0.9	± 0.3	
CP44.7*	2220-CP		44.7	± 0.9	± 0.3	
CP50	2222-CP		50	± 0.9	± 0.3	
CP60*	2224-CP		60	± 0.9	± 0.3	
CP70*	2226-CP		70	± 0.9	± 0.3	
CP75*	2228-CP		75	± 0.9	± 0.3	
CP80*	2230-CP		80	± 0.9	± 0.3	
CP90*	2232-CP		90	± 0.8	± 0.3	
CP100	2234-CP		100	± 0.8	± 0.3	
CP120*	2236-CP		120	± 0.8	± 0.3	
CP150*	2238-CP	150	± 0.7	± 0.3		
CP200	2240-CP	200	± 0.6	± 0.2	1000 μL	
CP200A*	2242-CP	200	± 0.6	± 0.2		
CP220*	2244-CP	220	± 0.6	± 0.2		
CP250	2246-CP	250	± 0.6	± 0.2		
CP300	2248-CP	300	± 0.6	± 0.2		
CP400*	2250-CP	400	± 0.6	± 0.2		
CP450*	2252-CP	450	± 0.6	± 0.2		
CP500	2254-CP	500	± 0.6	± 0.2		
CP600*	2256-CP	600	± 0.6	± 0.2		
CP700*	2258-CP	700	± 0.6	± 0.2		
CP750*	2260-CP	750	± 0.6	± 0.2		
CP800*	2262-CP	800	± 0.6	± 0.2		
CP900*	2264-CP	900	± 0.6	± 0.2		
CP1000	2266-CP	1000	± 0.6	± 0.2		

*Special versions upon request