

BIOHIT pipettes			mLine single channel					mLine multichannel			Proline single channel						Proline multichannel			Proline Plus single channel					Proline Plus multichannel										
HTL tips	Tip volume [μl]	Cat number	0,1-3	0,5-10	2-20	10-100	20-200	100-1000	0,5-10	5-100	30-300	0,1-2,5	0,5-10	2-20	5-50	5-100	10-100	20-200	100-1000	0,5-10	5-50	50-300	0,1-3	0,5-10	2-20	10-100	20-200	100-1000	0,5-10	10-100	30-300				
NON FILTER	0,5-10	30011 34011	●	●					●			●	●								●			●	●						●				
		30012 34012	●	●									●									●			●	●						●			
		30013 34013	●	●									●												●	●						●			
		30014 34014	●	●									●										●			●	●					●			
		30015 34015	●	●									●										●			●	●					●			
	1-200	30201 30221					●	●			●	●		●	●	●	●	●	●				●	●		●	●	●					●		
		30202 30222					●	●					●	●		●	●	●	●				●			●	●	●					●		
		30203 30223					●	●					●	●		●	●	●	●				●			●	●	●					●		
		30224					●	●					●	●		●	●	●	●				●			●	●	●					●		
		30225					●	●					●	●		●	●	●	●				●			●	●	●					●		
	1-300	30301					●	●			●		●	●	●	●	●	●	●				●	●		●	●						●		
		30302					●	●			●		●	●		●	●	●	●				●	●		●	●						●		
		30303					●	●			●		●	●		●	●	●	●				●	●		●	●						●		
	100-1000	31001						●																								●			
		31002						●																								●			
		31003						●																								●			
	FILTER	0,5-10	30016 34016		●								●	●									●			●							●		
		0,5-20	30026			●								●																					
0,5-50		30056												●								●													
1-100		30106				●			●				●	●	●	●										●						●			
1-200		30206					●			●			●										●												
1-300		30306								●																									
100-1000		31006						●																								●			

● compatible

● incompatible