

SOCOREX pipettes			Accura Manual single channel							Accura Manual multichannel				Accura Manual XS single channel										
HTL tips	Tip volume [μl]	Cat number	0,1-2	0,5-10	1-10	2-20	5-50	10-100	20-200	100-1000	0,5-10	5-50	20-200	40-350	0,1-2	0,5-10	1-10	2-20	5-50	10-100	20-200	100-1000		
NON FILTER	0,5-10	30011 31011	●	●	●						●				●	●	●							
		30012 34012	●	●	●							●				●	●	●						
		30013 34013	●	●	●							●				●	●	●						
		30014 34014	●	●	●							●				●	●	●						
		30015 34015	●	●	●							●				●	●	●						
	1-200	30201 30221				●			●	●		●	●	●					●			●	●	
		30202 30222				●			●	●		●	●						●			●	●	
		30203 30223				●			●	●		●	●						●			●	●	
		30224				●			●	●		●	●						●			●	●	
		30225				●			●	●		●	●						●			●	●	
	1-300	30301							●	●		●	●	●							●	●	●	
		30302							●	●		●	●								●	●	●	
		30303							●	●		●	●								●	●	●	
	100-1000	31001									●													●
		31002									●													●
		31003									●													●
	FILTER	0,5-10	30016 34016	●	●	●						●				●	●	●						
		0,5-20	30026			●	●												●					
		0,5-50	30056					●				●								●				
		1-100	30106						●													●		
1-200		30206							●			●	●									●		
1-300		30306																						
100-1000		31006								●													●	

● compatible

● incompatible