

DISCOVERY COMFORT VARIABLE VOLUME PIPETTE SERIES

SINGLE CHANNEL PIPETTES							MULTICHANNEL PIPETTES								
MODEL	CAT. NO.	COLOUR	VOLUME (µL)	A (%)	P (%)	FIT TO TIPS	MODEL	CAT. NO.	COLOUR	VOLUME (µL)	A (%)	P (%)	FIT TO TIPS		
DV2* D2	4041-DV 4051-D	RED	0.2	± 12.0	± 6.0	10µL	DV8-10 DV12-10	5121 5125	RED	MIN 0.5	± 10.0	± 8.0	10µL		
			1.0	± 2.7	± 1.3					5	± 4.0	± 2.0			
			MAX 2.0	± 1.5	± 0.7					MAX 10	± 2.0	± 1.2			
DV10 D10	4042-DV 4052-D	RED	MIN 0.5	± 4.0	± 2.8		200µL	DV8-50 DV12-50	5122 5126	YELLOW	MIN 5	± 4.0	± 2.5	200µL	
			5	± 1.0	± 0.6						25	± 3.0	± 1.2		
			MAX 10	± 0.5	± 0.4						MAX 50	± 1.6	± 0.6		
DV20 D20	4043-DV 4053-D	YELLOW	MIN 2	± 3.0	± 1.5			DV8-200 DV12-200	5123 5127	YELLOW	MIN 20	± 3.0	± 1.5		300µL
			10	± 1.0	± 0.5						100	± 1.5	± 0.8		
			MAX 20	± 0.8	± 0.3						MAX 200	± 1.0	± 0.6		
DV50 D50	4047-DV 4057-D	YELLOW	MIN 5	± 2.5	± 2.0		DV8-300 DV12-300	5124 5128	GREEN	MIN 50	± 1.6	± 1.5	300µL		
			25	± 1.0	± 0.6	150				± 1.2	± 1.0				
			MAX 50	± 0.8	± 0.4	MAX 300				± 1.0	± 0.6				
DV100 D100	4044-DV 4054-D	YELLOW	MIN 10	± 1.6	± 0.8										
			50	± 0.8	± 0.24										
			MAX 100	± 0.8	± 0.2										
DV200 D200	4045-DV 4055-D	YELLOW	MIN 20	± 1.2	± 0.6										
			100	± 0.8	± 0.25										
			MAX 200	± 0.6	± 0.2										
DV250 D250	4040-DV 4050-D	GREEN	MIN 50	± 1.0	± 0.4										
			125	± 0.8	± 0.3										
			MAX 250	± 0.6	± 0.3										
DV1000 D1000	4046-DV 4056-D	BLUE	MIN 100	± 1.6	± 0.4										
			500	± 0.7	± 0.2										
			MAX 1000	± 0.6	± 0.15										
DV5000 D5000	4048-DV 4058-D	WHITE	MIN 500	± 1.2	± 0.5										
			2500	± 0.6	± 0.2										
			MAX 5000	± 0.5	± 0.15										
DV10000 D10000	4049-DV 4059-D	WHITE	MIN 1000	± 2.5	± 0.6										
			5000	± 0.8	± 0.3										
			MAX 10000	± 0.5	± 0.2										

*DV/D2 - minimum volume 0.1 µl