

HTL MINI MICROCENTRIFUGE



HTL Mini Microcentrifuge is designed and dedicated to a wide range of applications in medical, veterinary, biochemical and other laboratories where constituents of investigated materials are to be separated under the effect of centrifugal force. The product is simple and convenient to operate. It is equipped with a quick-release rotor system and electronic brake. Once samples are loaded and the lid is closed, the rotor rapidly accelerates to 6,000 rpm, making it an ideal centrifuge for bringing small droplets to the bottom of tubes for micro-filtrations or basic separations.

- Accommodates up to 8 microcentrifuge tubes or 4 PCR strip tubes
- Comes with adapters for smaller tubes such as 0.2, 0.25, and 0.5 mL
- Quick release rotor system – no tools needed to exchange rotors
- Electronic braking - allows for more rapid deceleration
- Dual lid switch
- Compact, rugged design
- 2-year warranty



Parameters:

- Maximum Capacity: 8x1.5/2.0mL or 4x0.2mL strips
- Maximum speed: 6,000rpm
- RCF: 2,000xg
- Dimensions (LxWxH): 5.9x5.9x4.6in. (15x15x11.6cm)
- Weight: 2.65 lbs./1.2kg
- Electrical: Universal 100-240V
- Compliance: CE/TUV/CSA

Accessories included in the box:

Product	Cat. no.
HTL Mini Microcentrifuge	8000
Rotor for 8 x 1.5/2.0mL tubes	6770-RT
Rotor for 4 x 0.2mL strip tubes	6770-RTS
Individual adapters for 0.2mL tubes, fits in rotor 6770-RT, pack of 8	6770-RTA.20
Individual adapters for 0.25mL tubes, fits in rotor 6770-RT, pack of 8	6770-RTA.25
Individual adapters for 0.5mL tubes, fits in rotor 6770-RT, pack of 8	6770-RTA.5
Power adapter, 24V	401025



PZ HTL S.A.

Daniszewska 4 ■ 03-230 Warsaw ■ POLAND

Tel.: +48 22 492 19 00 ■ Fax: +48 22 492 19 93 ■ htlcs@corning.com

www.htl.com.pl