

PIPETY BRAND			Transferpette jednokanałowa										Transferpette wielokanałowa							Transferpette S jednokanałowa							Transferpette S wielokanałowa				
KOŃCÓWKI HTL	Zakres pojemności końcówek [μl]	Nr katalogowy	0,1-1	0,5-10	2-20	5-50	10-100	20-200	25-250	100-1000	500-5000	0,5-10	2-20	2,5-25	5-50	10-100	20-200	30-300	0,1-1	0,5-10	2-20	10-100	20-200	30-300	100-1000	0,5-10	5-50	10-100	20-200	30-300	
			BEZ FILTRA	0,5-10	301-03-052, 301-22-052	●	●								●							●							●		
301-03-102, 301-22-102	●	●											●							●							●				
301-03-152, 301-22-152	●	●											●							●							●				
301-03-351	●	●											●							●							●				
301-79-102	●	●											●							●							●				
1-200	301-02-052, 301-06-052				●	●	●	●	●				●	●	●	●	●	●			●	●	●				●	●	●	●	●
	301-02-152				●	●	●	●	●				●	●	●	●	●	●			●	●	●				●	●	●	●	●
	301-02-202, 301-06-352				●	●	●	●	●				●	●	●	●	●	●			●	●	●				●	●	●	●	●
	301-06-252				●	●	●	●	●				●	●	●	●	●	●			●	●	●				●	●	●	●	●
	301-79-052				●	●	●	●	●				●	●	●	●	●	●			●	●	●				●	●	●	●	●
1-300	301-16-052						●	●	●						●	●	●				●	●						●	●	●	●
	301-16-102						●	●	●						●	●	●				●	●						●	●	●	●
100-1000	301-01-202										●														●						
	301-01-302										●														●						
	301-01-352									●														●							
500-5000	301-71-052									●																					
Z FILTREM	0,5-10	302-05-102, 302-26-152	●	●								●	●						●							●					
	0,5-20	302-03-102			●								●							●											
	1-100	302-08-102					●							●							●						●				
	1-200	302-04-102						●	●						●						●							●			
	1-300	302-12-102							●								●											●			
	100-1000	302-01-102								●								●												●	

- kompatybilne
- niekompatybilne