

PIPETY CAPP			CappAero Comfort jednokanałowa					CappAero wielokanałowa					ecopipette jednokanałowa					CappTronic jednokanałowa			CappTronic wielokanałowa								
KOŃCÓWKI HTL	Zakres pojemności końcówek [μl]	Nr katalogowy	0,5-10	2-20	5-50	10-100	20-200	100-1000	0,5-10	2-20	5-50	10-100	20-200	30-300	0,2-2	0,5-10	2-20	5-50	10-100	20-200	100-1000	2-20	10-200	100-1000	2-20	10-200	100-1000		
			BEZ FILTRA	0,5-10	301-03-052, 301-22-052	●												●	●										
301-03-102, 301-22-102	●								●							●	●												
301-03-152, 301-22-152	●								●							●	●												
301-03-351	●								●							●	●												
301-79-102	●								●							●	●												
1-200	301-02-052, 301-06-052					●	●	●				●	●					●	●	●			●			●		●	
	301-02-152					●	●	●				●	●	●				●	●	●			●			●		●	
	301-02-202, 301-06-352					●	●	●				●	●	●				●	●	●			●			●		●	
	301-06-252					●	●	●				●	●	●				●	●	●			●			●		●	
	301-79-052					●	●	●				●	●	●				●	●	●			●			●		●	
1-300	301-16-052					●	●					●	●	●					●	●			●			●		●	
	301-16-102					●	●					●	●	●					●	●			●			●		●	
100-1000	301-01-202								●													●			●			●	
	301-01-302								●													●			●			●	
	301-01-352								●													●			●			●	
Z FILTREM	0,5-10	302-05-102, 302-26-152	●												●														
	0,5-20	302-03-102																											
	1-100	302-08-102				●					●								●										
	1-200	302-04-102					●					●								●			●			●			
	1-300	302-12-102											●																
	100-1000	302-01-102						●													●			●			●		

● kompatybilne

● niekompatybilne