

BRAND pipettes			Transferpette single channel								Transferpette multichannel						Transferpette S single channel						Transferpette S multichannel									
HTL tips	Tip volume [µl]	Cat number	0,1-1	0,5-10	2-20	5-50	10-100	20-200	25-250	100-1000	500-5000	0,5-10	2-20	2,5-25	5-50	10-100	20-200	30-300	0,1-1	0,5-10	2-20	10-100	20-200	30-300	100-1000	0,5-10	5-50	10-100	20-200	30-300		
NON FILTER	0,5-10	301-03-052, 301-22-052	●	●								●								●							●					
		301-03-102, 301-22-102	●	●									●								●							●				
		301-03-152, 301-22-152	●	●									●								●							●				
		301-03-351	●	●									●								●							●				
		301-79-102	●	●									●								●							●				
	1-200	301-02-052, 301-06-052			●	●	●	●	●					●	●	●	●	●	●			●	●	●					●	●	●	●
		301-02-152			●	●	●	●	●					●	●	●	●	●	●			●	●	●					●	●	●	●
		301-02-202, 301-06-352			●	●	●	●	●					●	●	●	●	●	●			●	●	●					●	●	●	●
		301-06-252			●	●	●	●	●					●	●	●	●	●	●			●	●	●					●	●	●	●
		301-79-052			●	●	●	●	●					●	●	●	●	●	●			●	●	●					●	●	●	●
	1-300	301-16-052					●	●	●							●	●	●				●	●						●	●	●	●
		301-16-102					●	●	●							●	●	●				●	●						●	●	●	●
	100-1000	301-01-202									●															●						
		301-01-302									●															●						
		301-01-352									●															●						
	500-5000	301-71-052									●																					
	FILTER	0,5-10	302-05-102, 302-26-152	●	●								●	●							●							●				
0,5-20		302-03-102			●								●								●											
1-100		302-08-102					●								●							●						●				
1-200		302-04-102						●	●						●							●						●				
1-300		302-12-102							●						●													●			●	
100-1000		302-01-102								●																						

● compatible

● incompatible