

GILSON pipettes			Pipetman Classic single channel					Pipetman G single channel					Pipetman L single channel					Pipetman L multichannel				Pipetman Neo single channel					Pipetman Neo multichannel								
HTL tips	Tip volume [μl]	Cat number	0,2-2	1-10	2-20	20-100	50-200	200-1000	0,2-2	1-10	2-20	10-100	20-200	100-1000	0,2-2	1-10	2-20	10-100	20-200	100-1000	0,5-10	2-20	20-200	20-300	0,2-2	1-10	2-20	10-100	20-200	100-1000	2-20	20-200			
NON FILTER	0,5-10	301-03-052, 301-22-052	●	●					●	●					●	●									●	●									
		301-03-102, 301-22-102	●	●					●	●					●	●									●	●									
		301-03-152, 301-22-152	●	●					●	●					●	●									●	●									
		301-03-351	●	●					●	●					●	●									●	●									
		301-79-102	●	●					●	●					●	●									●	●									
	1-200	301-02-052, 301-06-052			●	●	●					●	●	●				●	●	●			●					●	●	●			●	●	
		301-02-152			●	●	●					●	●	●				●	●	●								●	●	●					
		301-02-202, 301-06-352			●	●	●					●	●	●				●	●	●								●	●	●					
		301-06-252			●	●	●					●	●	●				●	●	●								●	●	●					
		301-79-052			●	●	●					●	●	●				●	●	●								●	●	●					
	1-300	301-16-052				●	●						●	●					●	●			●					●	●	●			●	●	
		301-16-102				●	●						●	●					●	●									●	●				●	
	100-1000	301-01-202						●							●						●										●				
		301-01-302						●							●						●										●				
		301-01-352						●							●						●										●				
	FILTER	0,5-10	302-05-102, 302-26-152	●	●					●	●					●	●									●	●								
		0,5-20	302-03-102			●						●						●										●						●	
		1-100	302-08-102				●						●						●										●						
1-200		302-04-102					●						●						●			●							●				●	●	
1-300		302-12-102																																	
100-1000		302-01-102						●							●						●										●				

● compatible

● incompatible