
















VIEW	VOLUME OF TIP [ $\mu$ L]	PIPETTE VOLUME RANGE [ $\mu$ L]	CAT. NO.	STERILE	PACKAGING
<b>STANDARD TIPS</b>					
	0.1 - 10 SHORT	0.1 - 2.0 0.5 - 10	301-03-052	NO	1000 TIPS/BAG
			301-03-102	NO	96 TIPS/BOX x10 BOXES
			301-03-152	YES	96 TIPS/BOX x10
			301-79-102	NO	10 REMOVABLE TRAYS OF 96 TIPS*
	0.1 - 10 EXTENDED	0.1 - 2.0 0.5 - 10	301-22-052	NO	1000 TIPS/BAG
			301-22-102	NO	96 TIPS/BOX x10
			301-22-152	YES	96 TIPS/BOX x10
	1 - 200 WITHOUT GRADUATION MARKS	0.5 - 20 5 - 50 10 - 100 20 - 200	301-02-052	NO	1000 TIPS/BAG
			301-02-152	NO	96 TIPS/BOX x10
			301-02-202	YES	96 TIPS/BOX x10
	1 - 200 WITH GRADUATION MARKS AT 10, 50, 100 $\mu$ L	0.5 - 20 5 - 50 10 - 100 20 - 200	301-06-052	NO	1000 TIPS/BAG
			301-06-352	YES	96 TIPS/BOX x10
			301-79-052	NO	10 REMOVABLE TRAYS OF 96 TIPS*
	1 - 300 FINE-POINT	50 - 250 50 - 300	301-16-052	NO	1000 TIPS/BAG
			301-16-102	NO	96 TIPS/BOX x10
	100 - 1000 BEVELLED	100 - 1000	301-01-202	NO	1000 TIPS/BAG
			301-01-302	NO	100 TIPS/ BOX x10
			301-01-352	YES	100 TIPS/ BOX x10
	500 - 5000 MACROVOLUME	500 - 5000	301-71-052	NO	250 TIPS/BAG
	1000 - 10000 MACROVOLUME	1000 - 10000	301-75-052	NO	200 TIPS/BAG
* REFILL CONTAINS 10 TRAYS WITH TIPS PLUS 2 RACKS					
<b>FILTER TIPS</b>					
	0.1 - 10 SHORT	0.1 - 2.0 0.5 - 10	302-05-102	YES	96 TIPS/BOX x10
	0.1 - 10 EXTENDED	0.1 - 2.0 0.5 - 10	302-26-152	YES	96 TIPS/BOX x10
	1 - 20 GRADUATION MARKS	0.5 - 20	302-03-102	YES	96 TIPS/BOX x10
	1 - 100 GRADUATION MARKS	10 - 100	302-08-102	YES	96 TIPS/BOX x10
	1 - 200 GRADUATION MARKS	20 - 200	302-04-102	YES	96 TIPS/BOX x10
	1 - 300 FINE-POINT	50 - 250 50 - 300	302-12-102	YES	96 TIPS/ BOX x10
	100 - 1000 BEVELLED	100 - 1000	302-01-102	YES	100 TIPS/BOX x10

© 2021 Corning Incorporated. All Rights Reserved.

February, 2021

© 2021 Corning Incorporated. All Rights Reserved.