

Clinipet⁺

Fixed Volume Pipettors

Clinipet⁺ Fixed Volume Pipettors are available in 17 models ranging from 5 µL up to 1,000 µL. The handle and mechanism are specially designed for fixed volume pipettors. The pipettors offer reduced pipetting forces, excellent accuracy and precision and a highly durable, UV-resistant and fully autoclavable body.



- Recalibration system allows for quick and accurate pipettor calibration

fully
121°C
autoclavable

UV
resistant

years
3
warranty

- Low pipetting forces prevent RSI
- Highly durable PVDF shafts guarantee optimal mechanical and chemical resistance



- Unique tip sealing system



Sliding tip ejector button – ideal for both right- and left-handed users.

Unique design to reduce tip ejection forces.



Easy removable tip ejector enables easy work with narrow tubes.

Clinipet+ Fixed Volume Pipettors

Model	Cat. No.	Color Code	Volume (μL)	A (%)	P (%)	Non-filter Tips
CP5	2204-CP	●	5	±2.0	≤1.0	200 μL
CP10	2208-CP	●	10	±1.2	≤0.5	
CP20	2212-CP	●	20	±0.9	≤0.4	
CP25	2214-CP	●	25	±0.9	≤0.4	
CP30*	2216-CP	●	30	±0.9	≤0.3	
CP50	2222-CP	●	50	±0.9	≤0.3	
CP75*	2228-CP	●	75	±0.9	≤0.3	
CP100	2234-CP	●	100	±0.8	≤0.3	
CP150*	2238-CP	●	150	±0.7	≤0.3	
CP200	2240-CP	●	200	±0.6	≤0.2	
CP200A*	2242-CP	●	200	±0.6	≤0.2	1,000 μL
CP220*	2244-CP	●	220	±0.6	≤0.2	
CP250	2246-CP	●	250	±0.6	≤0.2	
CP300	2248-CP	●	300	±0.6	≤0.2	
CP500	2254-CP	●	500	±0.6	≤0.2	
CP750*	2260-CP	●	750	±0.6	≤0.2	
CP1000	2266-CP	●	1,000	±0.6	≤0.2	

*Special version upon request.



Plexi 4-position stand
(Cat. No. 5480)

Carousel stand
(Cat. No. 5484)

Multiple stand
(Cat. No. 5449)

Shelf clip
(Cat. No. SP19451)