

Pipette Tips

Standard and filtered

Pipette tips are integral components of any pipetting system. Their shape, material properties, and purity have an influence on the pipetting results. These pipette tips are manufactured from virgin polypropylene in clean room conditions and free from biocontaminants. The range includes standard, sterile, and filtered versions in all standard packaging formats.

- ▶ For use with HTL™ pipettors for optimal accuracy and precision
- ▶ Compatible with many pipettor brands
- ▶ Strict quality standards (pursuant to ISO 9001) followed from R&D to production and delivery
- ▶ Premium quality virgin polypropylene used as raw material
- ▶ Suitable for sterilization with autoclave (121°C, 15 minutes, without dry cycle)

Bulk Packs

Pipette tips are available in bulk packaging, resealable plastic bags – keeping them safe from contamination. One bag contains 200, 250 or 1,000 tips depending on the tip volume.



Racks

Racked pipette tips are in a durable polypropylene box for good stability on the lab bench. One rack contains 96 or 100 tips depending on the tip volume.



Refill System

One refill system contains 10 trays with 96 tips plus an empty rack box. Available only for 10 µL and 200 µL tips with graduation marks.



Standard pipette tips

Cat. No.	Volume (µL)	Pipettor Volume Range (µL)	Sterile	Packaging	Qty
301-03-052			No	Bulk	1,000 pipette tips/pack
301-03-102			No	Rack	96 pipette tips/box, 10 boxes
301-03-152	0.1 - 10 (short)	0.1 - 2 0.5 - 10	Yes	Rack	96 pipette tips/box, 10 boxes
301-79-102			No	Refill	96 pipette tips/tray, 10 removable trays*
301-22-052			No	Bulk	1,000 pipette tips/pack
301-22-102	0.1 - 10 (extended)	0.1 - 2 0.5 - 10	No	Rack	96 pipette tips/box, 10 boxes
301-22-152			Yes	Rack	96 pipette tips/box, 10 boxes
301-02-052			No	Bulk	1,000 pipette tips/pack
301-02-152	1 - 200	2 - 20 5 - 50 10 - 100 20 - 200	No	Rack	96 pipette tips/box, 10 boxes
301-02-202			Yes	Rack	96 pipette tips/box, 10 boxes
301-06-052			No	Bulk	1,000 pipette tips/pack
301-06-352	1 - 200 (graduation marks at 10, 50, 100 µL)	2 - 20 5 - 50 10 - 100 20 - 200	Yes	Rack	96 pipette tips/box, 10 boxes
301-79-052			No	Refill	10 removable trays of 96 pipette tips*
301-16-052			No	Bulk	1,000 pipette tips/pack
301-16-102	1 - 300 (fine-point)	50 - 250 50 - 300	No	Rack	96 pipette tips/box, 10 boxes
301-01-202			No	Bulk	1,000 pipette tips/pack
301-01-302	100 - 1,000 (beveled)	100 - 1,000	No	Rack	100 pipette tips/box, 10 boxes
301-01-352			Yes	Rack	100 pipette tips/box, 10 boxes
301-71-052	500 - 5,000	500 - 5,000	No	Bulk	250 pipette tips/pack, 10 packs
301-75-052	1,000 - 10,000	1,000 - 10,000	No	Bulk	100 pipette tips/pack, 20 packs

*Refill contains 10 trays with pipette tips plus 2 racks.

Filtered pipette tips, Rack

Cat. No.	Volume (µL)	Pipettor Volume Range (µL)	Sterile	Qty
302-05-102	0.1 - 10 (short)	0.1 - 2 0.5 - 10	Yes	96 pipette tips/box, 10 boxes
302-26-152	0.1 - 10 (extended)	0.1 - 2 0.5 - 10	Yes	96 pipette tips/box, 10 boxes
302-03-102	0.1 - 20 (graduation marks)	2 - 20	Yes	96 pipette tips/box, 10 boxes
302-08-102	1 - 100 (graduation marks)	10 - 100	Yes	96 pipette tips/box, 10 boxes
302-04-102	1 - 200 (graduation marks)	20 - 200	Yes	96 tips/box, 10 boxes
302-12-102	1 - 300 (fine-point)	50 - 250 50 - 300	Yes	96 tips/box, 10 boxes
302-01-102	100 - 1,000 (beveled)	100 - 1,000	Yes	100 pipette tips/box, 10 boxes

