

HTL™ Swiftpet™ Pro Controller

Pipet Controller

Pipetting using glass and plastic volumetric pipets has never been easier than with the Swiftpet Pro Pipet Controller. An ergonomically shaped handle along with optimally located, smooth pushbuttons and switches offer effortless pipetting even during periods of extended use. A convenient charging stand allows the pipet controller to be stored when not in use. Easily accessible switches on the Swiftpet Pro allow the user to select different operation modes depending on the volume of pipet and viscosity of liquid. The Swiftpet Pro is an excellent tool for pipetting in tissue culture laboratories, water treatment facilities, academic labs, and more.

Safe and efficient work

- ▶ Protected by a PTFE filter that blocks liquids from entering the unit
- ▶ Autoclavable filter, pipet holder and nose piece
- ▶ UV-resistant body for safe sterilization
- ▶ Compliance: CE/UKCA



Comfort while working

- ▶ An ergonomically shaped handle
- ▶ A convenient wall hanger allows the pipet controller to be stored when not in use
- ▶ Suitable for glass and plastic volumetric pipets 0.5 - 100 mL
- ▶ Powerful Li-ion battery enable 6 hours of continuous work



- ▶ LCD display showing battery charge level



Adjustment of speed and working mode

- ▶ Easily accessible function buttons for control of speed and working mode
- ▶ Additional speed adjustment via the pressure applied to the pushbuttons
- ▶ LCD display showing battery charge level

Product	Cat. No.
Swiftpet Pro	4743

Accessories and Spare Parts	Cat. No.
PTFE autoclavable filter, 0.2 µm, pack of 5	SP9143
PTFE autoclavable filter, 0.45 µm, pack of 5	SP9144
Charging stand, 2-positions	SP29056
Charging stand, 1-position	5411
Wall hanger	SP9029
Silicone pipet holder	SP29054
Nose piece	SP29053
Li-ion battery set	SP29096
Battery cover	SP29090
Universal power supply, 9V with set of adapters (EU, USA, GBR, KOR, and AUS)	SP29100